



# OVERVIEW



# INSTALLATION



## Fasten Sensor Unit

First, remove the mattress from the bed. Then select the immovable bed part, which is usually situated in the middle of the bed.

Attach the sensor unit on the immovable bed part, using the straps attached to the sensor unit.

**Important:** Make sure the labelling points upwards and that the sensor unit cannot be moved.



## Check Cable Routing

Attach the sensor cable to the bed frame. The supplied cable ties can be used. Make sure that the sensor cable can't get stuck between bed parts.

**Important:** Move bed parts and check if nothing gets blocked!



## Test Nurse Call

Attach the handheld device to the bed and connect the sensor unit. Plug in the power supply. Connect the nurse call cable (optional).

**Important:** Please check if the nurse call is working. Press the "Keylock" (🔒) and then press "Bed Exit Alert" and "Mobility Alert" (🚰 + 🚰) at once for 2s.



## Calibration / Tune to Mattress

The sensor unit must be fine-tuned to the mattress in order to detect movements correctly. Unlock the keys (1.) and calibrate the device (2.).

1. Press "Keylock" (🔒)
2. Press "Keylock" and "Repositioning Button" (🔒 + 🔄) at once for 4s, until the traffic light flashes and the calibration process starts.

**Important:** The mattress must be placed on top of the sensor and the bed must not be occupied.

The calibration must be done after each installation and on every change of patient. The device must also be calibrated in case the sensor unit or handheld device is being replaced.



**Calibration successful**



**Calibration failed**