

Our mutual goal:
Satisfied patients and residents thanks to needs-oriented care


- ✓ Assessment
- ✓ Prevention of Falls & Pressure Ulcers
- ✓ Sleep Patterns

**Assure Quality &
Reduce Effort**



compliant concept AG
Undermülstrasse 28
CH-8320 Fehraltorf
Tel: +41 44 552 15 00
Fax: +41 44 552 15 09
info@compliant-concept.ch
www.compliant-concept.ch

mm mobility
monitor

 designed and made in Switzerland

The Modular Nursing Care Solution



The Mobility Monitor is an assessment instrument for nursing; care measures can be carried out in a more targeted way.




The Mobility Monitor functions are available as modules. The modules can be combined at will. Upgrades are possible at all times.

Assure Quality & Reduce Effort

- ✓ Assessment
- ✓ Prevention of Falls & Pressure Ulcers
- ✓ Sleep Patterns

The use of the Mobility Monitor is easy, noiseless and discreet. The measurement unit is placed under the mattress and connected to the control unit.



The lights on the console display the mobility of the patient in an easy way.



Efficient prevention of falls thanks to bed edge and bed exit warnings.



Targeted repositioning thanks to immobility warning when patient fails to move sufficiently.



The Mobility Monitor can be connected to all common nurse call systems.

Efficient Assessment & Purposeful Planning




Mobility data facilitates cooperation and fact-based discussion between medical and care staff.



Easy data transfer via USB stick or radio transmission.



Intuitive mobility profile displayed on screen.



Example of mobility analysis: Loss of mobility (01:50am) after taking medication. A nurse call triggered by the system (03:50am) allowed urgently necessary repositioning (white bar).